

# IONIZED OXYGEN INFUSER

CONSUMER  
BEAUTY  
TOOL

## Create Oxygenated Mist for Face and Body



**SOUND BITE:** *"Help skin breathe, from head to toe."*

**Problem:** Every day, oxygen levels in skin decline, bombarded by environmental stressors that trigger dehydration, wrinkles, sagging, uneven skin tone, blemishes, and inflammation.

**Solution:** IONIZED OXYGEN INFUSER dispenses oxygenated skincare nutrients in an ultra-fine mist that instantly hydrates, nourishes, and brightens skin on face and body. Oxygenated, micronized mist propels skincare product molecules into skin for an immediate dewiness and healthy glow. Easy-to-use dispenser leaves skin cool, refreshed, and rejuvenated.

**Use with water-based products only.**

### WHAT DOES IT DO?

- Saturates skin in oxygenated skincare mist
- Instantly hydrates & creates a dewy glow
- Cools, refreshes & rejuvenates face AND body
- Helps ward off effects of environmental stressors

### DIRECTIONS

Pour Oh My Glow Activator into mini-flask attached to top of infuser. Add Oh My Glow Serum. Close lid and gently shake flask to blend liquid. Press power button and apply oxygenated mist. Press power button to turn off machine.

### SKIN TYPES

all skin types

### CONDITIONS

Dryness, lines & wrinkles, uneven tone

### CAUTION

For use with water-based products only. Do not use when unit is plugged in. Do not use while unit is charging. Do not leave product inside unit for a prolonged period of time. Operate only according to intended use as described in user manual. Keep out of reach of children.