

BIOGRO LED BRUSH HAIR SERUM APPLICATOR

Infuse & Energize



SOUND BITE: "Give thinning hair the brush-off"

BioGro LED Brush is a multi-function, easy-to-use hair serum applicator brush that evenly distributes serum onto scalp to facilitate delivery of

Kizo Lab Follicle Activator Hair Serum.

Follicle-energizing tool features LED red light, LED purple light, and 2-speed massage to stimulate scalp and hair follicles, addressing hair loss and thinning for thicker, fuller-looking hair.

WHAT DOES IT DO?

- Evenly distributes hair serum onto scalp
- Energizes scalp and hair follicles
- Addresses hair loss and thinning
- Anti-inflammatory, scalp-stimulating red and purple LED

QUICK START

STEP 1: Fill cartridge with **Kizo Lab Follicle Activator Hair Serum**. Press and hold Power button for 3-4 seconds. Quickly press Mode button twice.

STEP 2: Starting at hairline, press brush against scalp so teeth are in contact with scalp. Glide brush across scalp from front to back until targeted areas have been covered. Use at least 3 times per week.

DETAILED INSTRUCTIONS

For detailed instructions, consult User Guide or lemieuxpro.com

CONDITIONS

Hair loss, thinning, brittleness, dullness, breakage

CAUTION

Keep out of reach of children. Operate this unit only as described in User Guide. This unit should not be used by pregnant women or individuals with heart disease, abnormal blood pressure, or cerebrovascular disease.