

## MASK

# Boto-Needle Frown Line Patch

Smooth & Plump

**sound bite:** "Precision infusion to target deep lines."

Bio-available micro needles dissolve into skin to reduce the appearance of forehead wrinkles and expression lines.

Fast-tracks the delivery of skin-firming peptides and line-plumping moisture-boosters, while stimulating skin and encouraging collagen synthesis for a more youthful look.

### what does it do?

Visibly targets forehead lines & wrinkles

Hydrates, smooths & plumps

Helps encourage collagen synthesis

Boosts absorption of underlying serum

### directions

Remove adhesive strip from edge of patch. Press patch gently onto cleansed, dry skin. For optimal results, apply serum and let dry before applying patch. (Note: if skin is wet before application, adhesive will not stick.) Leave patch on for a minimum of 30 minutes, or overnight.

### key ingredients

<b>Hyaluronic Acid</b>	Attracts & retains moisture to plump & improve barrier function
<b>Acetyl Hexapeptide-8</b>	Reduces appearance of wrinkles caused by repeated facial muscle movements
<b>sh-Oligopeptide-1</b>	Bio-identical, lab-synthesized peptide mimics EGF to promote skin repair & renewal
<b>Glutathione</b>	"Master Antioxidant" helps protect skin from signs of oxidative stress & damage

### ingredients

Oryza Sativa Bran Water, Hyaluronic Acid, Water, sh-Polypeptide-1, Niacinamide, Sodium Hyaluronate, Sodium Hyaluronate Crosspolymer, Glutathione, Ceramide EOP, Ceramide NS, Ceramide NP, Phospholipids, sh-Oligopeptide-1, sh-Oligopeptide-2, Acetyl Hexapeptide-8, Caprylyl Glycol, Ethylhexylglycerin, Undecylenoyl Glycine, Capryloyl Glycine.



### available in

→ 1 mask

### skin types

- Oily
- Balanced to Oily
- Balanced
- Balanced to Dry
- Dry

### skin concerns

- Fine lines
- Deep lines
- Loss of elasticity

### pH

N/A